



**Actual size: 140mm (Width)  
x 175mm (Height)**

Ruler is for scaling purposes.

**Use this fitting template as  
a guide to size your patients.**


**How to use the fitting template:**

- Cut out the template along the dotted line
- Follow directions on template to determine the recommended cushion size.
- Please ensure the template is printed to scale.

[www.EU-PAP.co.uk](http://www.EU-PAP.co.uk)




**Rx Only** 0123



Quattro™ Air  
AirFit™ F10

Full face mask  
Full face mask **for Her**

**L**  
**M**  
**S**  
**XS**  
**SIZE**



**Sample size check**

AirFit and Quattro are trademarks of ResMed Ltd.  
Quattro is registered in the U.S. Patent and Trademark Office.  
© 2014 ResMed Ltd. All rights reserved. Specifications may change without notice.  
Product availability may vary across regions.

**Rx Only** 0123

